



CMI

Calico Medical
INCORPORATED



X-ray Equipment & Accessories
Proudly Canadian 🇨🇦

Rayence Digital Flat Panels

- No hardware integration
- Uses existing cassette tray
- Auto trigger software
- Minimal installation time
- Available in both Gaddox (SGC) and Cesium (SCC) Scintillators



Calico Medical Inc. is proud to introduce the Rayence line of Digital Flat Panels to our product line.

This flat panel features 3.9lp/mm and a pixel pitch of 129 microns. The thin profile makes it an easy upgrade to any existing x-ray system.

If you are considering digital in the near future, it will be fast and seamless using the popular DX-R X-ray acquisition Software.

Setting up a demo in your office is a phone call away. Call today to discuss what we can do for you!

**The
1717**



**Rayence
SGC/SCC
detectors**

Robust, Comprehensive, Customizable

This software is the core of any digital X-ray unit:

- Image processing can be adapted to Individual user needs
- Developed to enable organ specific optimization (guaranteeing the highest quality X-ray images)
- Integrated functionality simplifies your daily routine (i.e. radiographic positioning guide)
- Easily integrated with existing patient management systems (Worklist, DICOM Send, etc.)
- Diagnose X-Ray images within the acquisition software (integrated full viewer)
- Can be applied as a complete diagnostic workstation with the option to upgrade to a PACS (Picture Archiving and Communications System)
- Fully modular with several options that can be added as your practice grows
- Outstanding image quality

Veterinary Digital Imaging Solutions- Intra-Oral Sensor

Raynence- EzSensor P Intraoral Dental Sensor HDI 2000

EzSensor P is the industry leading intra-oral sensor, offering superior diagnostic images that are essential to successful dental treatments. Features unparalleled durability and image quality, EzSensor P delivers speed and efficiency while keeping patients comfortable.

Features:

- Clear, High Quality Images
- Ergonomic design: Slim and rounded corner exterior
- Time-Saving workflow
- Rigorous testing ensures durability and reliability



The EzRay, VATECH's intra-oral X-ray imaging system



Reliable Design for Easy & Stable Use

The ergonomically designed X-ray tube head can be positioned at any anatomical imaging position. In addition, the strong and stable arm produces distortion-free images and allows simple positioning.

Features:

- Sophisticated Design
- Creates a modern and comfortable atmospheric in your clinic
- Brings relief to patients by reducing anxiety and fear of discomfort
- Intuitive & easy to use operating panel
- Straightforward process for image capture with Intuitive Icons
- Detachable Control Panel
- Round Shape and Lightweight Tube head
- Easy Grip with Ergonomic Style

My Vet DR Equine Portable Package

Dr Equine Portable Package

The Rayence DR portable package provides the ultimate in adaptability to its environments. It uses auto-triggered signal sensing technology that does not require any integration into an X-ray generator. The package comes with the Rayence XMARU 1012WCA Panel Dector, which is a digital detector that completely, wireless. For the convenient use of this unique wireless solution, size, thickness, and handle of this model is ergonomically designed. For the first time in the history of the industry that a function to install and remove the detector handle horizontally and vertically has come to light. The panel has a charging stand that allows the unit to be fully changed when need of use.



Features:

- Shock Resistant
- Cloud
- Wireless
- Auto Image Rotation
- Sleep mode
- Light Weight
- Built in Image Storage
- Water Resistant
- Superior Image Quality
- Charging Stand



VET-20BT Hybrid Battery Portable

The VET-20BT Hybrid battery Portable is the lightest, most compact hybrid battery powered generator and the most advanced portable X-ray unit in the radiology market today. Vet-20BT hybrid battery X-ray unit is integrated with Wireless sync. module (optional for DR) which is facilitated by Radio-Frequency (RF) communication device.

Features:

- Light Weight
- Portable
- Lithium Ion Polymer Battery

MyVet Table for Small Animals

The MyVet's vet table is a revolutionary Veterinary Elevating X-Ray Table, which delivers silent elevating to fulfill comfort and stability. The MyVet's vet table is a cost-effective, workflow-efficient approach to digital radiographic imaging, designed to improve patient workflow in veterinary practice of small animals.



Sedecal Vet Ray X-Ray System

Sedecal is committed to developing the ultimate X-ray. Manufactured by Sedecal, systems can be customized.

Systems are compatible for Companion, Equine and mixed animal practices.

Conventional units are available for varying budgets.

Raymax DR Deluxe Model

The Raymax Veterinarian System features digital radiography (DR), a high frequency generator, and a 4 way float table.

It is an affordable complete Digital Radiography (DR) solution with superior quality.

Available in non-digital system as well.



SiMPLi Series USB Ultrasound

The SiMPLi Series ultrasound is the most affordable ultrasound unit that can be purchased today. It is a lightweight unit that plugs right into a USB port and is fully portable. You are able to plug it into any windows computer, laptop or tablet. The unit provides instant real time imaging right after you start.

Features:

- Lightweight
- Portable
- Real Time Imaging
- USB Port (any laptop or computer)
- Handheld
- Two sizes (142mm x 58mm and 150mm x 62mm)



Laptop not included

E-Cube i7 Portable Ultrasound

The E-CUBE i7 Veterinary ultrasound system is a cart-based system which has a powerful and qualified system architecture and user-optimized ergonomically design with a sliding keyboard cover. It delivers remarkably excellent imaging performance from routine examinations to complex cardiac examinations at a veterinarian's office.

Features:

- Lightweight
- 1-hour battery life
- Compact Laptop style with a carrying case Speed
- Fast boot-up with integrated SSD(Less than 45 seconds)
- Scanning Comfort
- 15.6" width & height with a resolution(1600x900) LED monitor with IPS technology (Wider angle, better expression of colour/contrast, less eye strain)
- Solid Aluminum Shell Library Quiet
- Silent fan noise of 31dp (comparing to a library noise level of 30dp)





FRONTAL APRON SIZES			
SIZE	Actual Apron Width x Length	Shoulder Width	Max. Chest/ Hip
X-Small	20" X 34"	14"	34"
Small	22" X 35"	14½"	38"
Medium	24" X 37"	16"	42"
Large	24" X 39"	18½"	44"
X-Large	26" X 41"	20½"	46"
XX-Large	28" X 43"	21½"	50"

Code	Description
(Apron)	
VA1	Veterinary Frontal with Velcro
VA2	Veterinary Frontal with Belt & Buckle
VA3	Veterinary Frontal with Tie-Straps

Code	Description
(Thyroid Shield)	
TSR-IT	Thyroid Shield Sewn-on (Integrated)
TSRL	Standard Thyroid Shield Loose
TSRT	Standard Thyroid Shield Tethered
TSRB	Standard Thyroid Shield Tethered with Buckles

VA1 ADJUSTABLE FIT FRONTAL Velcro Closure

- Easy on and off convenience.
- Adjustable Velcro arms provide excellent weight distribution and comfort.
- Adjustable shoulder strap prevents the apron from slipping.
- Optional no-charge chest pocket.

VA2 ADJUSTABLE FIT FRONTAL Belt and Buckle

- Easy on and off convenience.
- Adjustable heavy-duty 2" web belt with quick-release buckle closure.
- Adjustable shoulder strap prevents the apron from slipping.
- Optional no-charge chest pocket.

VA3 Economy Frontal Tie-Strap Closure

- Easy on and off convenience
- 1" straps cross the back and tie securely at the front
- Optional no charge chest pocket



PET TRANSPORTER

- A safe and secure way to move patients easily.
- Sturdy fiberglass poles with rubber handles.
- Stain-resistant nylon with 2" webbing and adjustable buckle reinforcement.
- Code: VPT

VETERINARY SANDBAGS *with or without handles*

- Sandbags are rectangular in shape. An inner bladder is filled with sand and then inserted into a blue nylon or vinyl bag.



Code (Nylon)	Code (Vinyl)	Description
VS85-N	VS85-V	VET Sandbag - 5lb - 7" x 8"
VS87-N	VS87-V	VET Sandbag - 7lb - 7" x 12"
VS812-N	VS812-V	VET Sandbag - 12lb - 7" x 16"
VS85H-N	VS85H-V	VET Sandbag - 5lb - 7" x 8" with Handles
VS87H-N	VS87H-V	VET Sandbag - 7lb - 7" x 12" with Handles
VS812H-N	VS812H-V	VET Sandbag - 12lb - 7" x 16" with Handles
VS812S-N	VS812S-V	VET Sandbag - Snake - 12lb - 5" x 27"



#TA-1 ea. set
R & L - 5/8"
2 Initials - 3/8" or 3 Initials - 1/4"
Without Initials, order:
#TA-4

Markers are available in custom colors.
Visit our website for more details

Open Hand Shield

This shield resembles a paddle with a cuff. It is used to protect the top of the hand. One-size fits all with a lead-equivalency of 0.5mm.

CODE: GL-OH



Open Palm Mitts

These mitts are primarily designed for veterinary applications where fingers and thumbs are available to be used in a procedure. One-size fits all with a lead-equivalency of 0.5mm.

CODE: GL-OP



Open Palm with Thumb

Identical to the Open-Palm Mitt with a closed thumb. One-size fits all with a lead-equivalency of 0.5mm.

CODE: GL-OPT



Ducky Mitts

These mitts are designed for veterinary applications where fingers and thumbs must be protected. One-size fits all with a lead-equivalency of 0.5mm.

CODE: GL-DM



Closed Mitts with Thumb

Similar to the Ducky Mitt without the finger pouch. One-size fits all with a lead-equivalency of 0.5mm.

CODE: GL-DMT



* Different colours shown are for deliniation purposes only.

All glasses are 0.75mm L.E. Side-Shields (If equipped) are 0.35mm L.E. All glasses include a protective case and strap.
All glasses are available in Rx: Single, Bifocal and Progressive.



UltraLite

Stylish and comfortable wraparound frame.
High-impact resistant plastic for safety.
Colours: Black, Bronze, Carbon Fibre, Pink, Purple/Silver, Red, Silver and Yellow/Blue.



Maxi

(Mini, Regular or Large)

Unisex shape suits any face and gender.
Padded temple bars and nose bridge for extra comfort.
The most durable product in our line!
Colours: Black, Blue, Red and Silver.



Sport Wrap

Solid, sturdy and extremely comfortable wraparound glasses.
Built-in side-shields for safety.
Colours: Black/Yellow and Brown/Orange.



Basic 52 Wrap

Excellent peripheral vision with leaded glass side-shields.
Economical and lightweight with a spring hinge.
Comes with a silicone, adjustable non-slip uni-fit bridge for increased comfort.
Colours: Tortoise.



Regular



Large


Fitover

(Regular or Large)

Designed to fitover existing prescription glasses.
Black rubberized nylon frame.
Lightweight and comfortable.
Colours: Black and Tortoise.

Veterinary Positioning Aids

- Various Wedges
- Head Holders
- Body Cradles
- U-Blocks
- Pre-set Vet Kits
- Custom Pieces Available




MXR Podoblock

In combination with our Podoholder or Podo-Prot. box this unique system guarantees maximum distance and no fingers on the X-rays.

The height from ground to tip off hoff in Podoblock:
90 mm, 120 mm, 170 mm

The Topblock is easily interchangeable (Angle diversity from 45° to 60°).



Scandroid

Multi-functional Scan Tool

The ideal height adjustable still to transport to the job site or for clinic use. Always ready for use. Under the seat, safely also during transport, you have storage for all necessary scanning tools.

Position One Height: 23cm

Scanning the foot and lower leg

Position Two Height: 40 cm

Scanning front knee and middle parts of leg

Position Three Height: 50cm

Scanning the stifle and higher parts of the leg



If extra height is required we have an optional gas spring with a height adjustment from 56cm to 81cm.



CMI
Calico Medical
INCORPORATED

Calico Medical Inc.
A leader in Comprehensive Imaging
Sales and Service

- CMI is a comprehensive diagnostic imaging company.
- We provide professional sales support, skilled installation, relocations, repairs and ongoing maintenance to hospitals, veterinarians, clinics, and chiropractors throughout Canada
- Our team stands behind every one of our imaging systems, ready to guide and support you from site planning to ongoing maintenance and everything in between
- Our company is qualified and capable of ensuring your specific imaging needed are met

For more information about the variety of products that we have to offer, please visit our website at www.calicomedical.com. For more up to date information about offers and new products like us on Facebook (www.facebook.com/CMI-Calico-Medical-286497931783068) and follow our Instagram (@calicomedical).

