



Veterinarian Digital Radiography (DR) X-Ray

CSI Cesium Iodide

The Affordable DR Solution with Superior Quality



Let Raymax imaging equipment show you what lies beyond...

Raymax equipment meets all North American Regulatory Standards.



The Raymax Veterinarian System features digital radiography (DR), a high frequency generator, and a 4 way float table. It is an affordable complete Digital Radiography (DR) solution with superior quality.



Features / Specifications

Height: 31 inches

Length: 60 inches

Width: 28 inches

Cross Table Travel: 8 inches


Longitudinal Table Travel: 14 inches

Space Requirement: 6 x 9 feet

Vertical SID Fixed 40 inches

Flat Plate Detector:

- Fixed configuration cassette detector
- True 43 cm x 43 cm (17" x 17") image
- High resolution 100 μ m pixel pitch
- Up to 65,536 grey levels (16-bit ADC)
- Automatic Exposure Detection (AED)
- Direct deposition CsI scintillator (XRpad 4343 F)



Computer Workstation (minimal req.)

- 5.0 GHz Intel CPU
- 4 GB High Speed RAM
- 1 TB Hard Disk Drive @ 7200 RPM
- Independent PCI Video Graphics Card (1 GB)
- DVD/CD Drive @ 16 RW with Burner
- Windows 7 Operating System

Monitor Specifications (minimal req.)

- 2MB High Contrast 1920 x 1080
- 22.5 Inch Display - LED
- Touch Screen (Optional)

*Appearance and specifications subject to change without notice

DISTRUBUTED BY:



1789 Notre Dame Drive
St Agatha, ON N0B 2L0
519-885-3254
office@calicomedical.com
www.calicomedical.com