


# MAKE YOUR PRACTICE EXTRAORDINARY

Veterinary Intraoral X-Ray Imaging System



## **EzRay**<sup>™</sup> Wall Mount X-Ray *A Veterinary Dental System*

- Easy-To-Move Anti-Drift Arms
- Detachable Control Panel
- Easy-To-Read Panel Icons
- Lightweight Tube Head
- Ergonomic, Easy-Grip Handle
- High Frequency DC X-ray for Superior Digital Images



# THE EZRAY INTRAORAL X-RAY IMAGING SYSTEM CREATES A MORE SOPHISTICATED CLINICAL ATMOSPHERE



## SPECIFICATIONS


|                      |                                                                  |
|----------------------|------------------------------------------------------------------|
| Focal Spot           | 0.4mm                                                            |
| Tube Voltage         | Selectable 60 - 70kV                                             |
| Current              | 4 - 7mA                                                          |
| Exposure Time        | Selectable 0.04 - 2.0 sec                                        |
| Arm Length (Options) | 62" / 68" / 80"                                                  |
| Weight               | 56.9lb (25.8kg) / 58.0lb(26.3kg) / 61.1lb(27.7kg)                |
| Line Voltage         | 100 - 230V $\pm$ 10%                                             |
| Accessories          | Exposure hand switch, Cone:Round 8"/16" Double stud wall bracket |
| Option               | Cone: Round (12") / Rectangular (8" / 12")                       |

\* Specifications subject to change without prior notice \*

## RELIABLE DESIGN FOR EASY & STABLE USE


The sophisticated, ergonomically designed X-ray tube head mounted on a strong and stable arm can be easily and comfortably positioned to acquire motion-free intraoral images creating stress-free atmosphere for animals in your clinic.

- Intuitive and easy-to-use operating panel
- DICOM Compliant System



## OBTAIN HIGH QUALITY DIAGNOSTIC IMAGES

With a 0.4mm focal spot, the EzRay intraoral X-ray system provides high resolution images for more exact diagnoses when acquiring maximum image information with the smallest focal spot technology.





## DIMENSIONS and OPTIONS

Note 1.  
Horizontal Arm Length Options  
• 17.7" • 23.6" • 35.4"

Note 2  
• 17.7" Horizontal Arm - 28.1"  
• 23.6" Horizontal Arm - 34"  
• 35.4" Horizontal Arm - 45.8"

Note 3  
• 17.7" Horizontal Arm - 62.4"  
• 23.6" Horizontal Arm - 68.3"  
• 35.4" Horizontal Arm - 80.1"  
(with 8" Cone)



**CMI**  
*Calico Medical*  
INCORPORATED

1789 Notre Dame Drive, Unit 103  
St Agatha, ON N0B 2L0  
519-885-3254 [www.calicomedical.com](http://www.calicomedical.com)

**EzRay**<sup>TM</sup>

