


Advantage Plus

by **TXR**

CM-CH1

17" x 17" UPRIGHT SYSTEM



The CM upright system is ideally suited for film or CR imaging. The system was designed to deliver high quality performance at an affordable price. Design criteria included improved esthetics as compared to traditional floor to wall or floor to ceiling mount tube stands.

Extensive travel of the x-ray tube and receptor provides safe imaging of pediatrics without the use of a stool while being able to obtain full coverage on tall patients.

The CM-CH1 is a sound investment that will deliver years of trouble-free service.

STANDARD FEATURES INCLUDE

FLOOR MOUNT TUBE STAND

- Platform tube mount, 90° tube rotation
- Magnetic locks
- 86.6" floor track provides 78.7" longitudinal tube travel

WALL STAND

- Floor to wall mount
- Magnetic locks
- 17" x 17" Bucky
- 103 Line, 8:1 ratio grid
- Cassette tray
- 59" Vertical travel

COLLIMATOR

- Manual
- Laser light
- Swivel mount

OPTIONS

- Choice of 6 generators
- AEC (Automatic Exposure Control)
- Various grids
- Various tubes

