




CMI
Calico Medical
INCORPORATED



Chiropractic Direct Digital X-Ray



**Let Raymax imaging equipment
show you what lies beyond...**

Raymax equipment meets all North American Regulatory Standards



CMI
Calico Medical
INCORPORATED

- Operating Modes: kV/mAs, kV/mA/Time
- 32/40/50/65/80 & 100 kW power levels
- Compatible with both digital and conventional radiographic (RAD) systems
- kV/mAs Stepping
- Entire apparatus self-supporting
- Space requirement 5 x 10 feet
- No architectural changes required
- Fast, easy installation
- Operator station shielding incorporated



Generator:

Generator Type

CPI 40 kW
High Frequency
Output up to 400 kHz
40 - 150 kV
10 - 500 mA
1 Phase, 230 VAC
768 Techniques
kV/mAs

kV (Radiographic)

mA Range RAD

Line Nominal/Phase²

Anatomical Programs RAD

Technique Selection



Approvals:

CSA, U/L, FDA, Health Canada

DISTRIBUTED BY:



CMI
Calico Medical
INCORPORATED

Calico Medical Inc.
1789 Notre Dame Dr, Unit 103, St. Agatha, Ontario Canada N0B 2L0
Phone: (519) 885-3254 | Email: office@calicomedical.com